

HAND & FOOT

SURGERY INSTRUMENTATION



- Soft Tissue Instruments
- Bone Instruments
- Tendon Instruments
- Small Joint Arthroscopy Instruments
- Micro Instruments



International Instruments
Designers & Manufacturers of Surgical Instruments



International Instruments

Designers & Manufacturers of Surgical Instruments

US FDA, CE & ISO REGD.

ABOUT US

JJ International Instruments - A leader in designing and manufacturing of high quality stainless steel surgical instruments in India. With a vast experience of 46 years in the industry, JJ provides unmatched service and unbeaten quality to the worldwide market. Corporate office and manufacturing unit with state of art facilities is located in south Indian state, Kerala - The renowned 'God's Own Country'.

CONTINUES SUCCESS...!

JJ became world wide since exhibited in AAOS 2013 with tremendous response from both domestic and international surgeons. JJ now participates many tradeshow in USA and more than 30 tradeshow all over India yearly. JJ is the only one company now in India who has approved by US FDA.

*" International
Quality
@
Indian Price "*

CE Reg No. DCS/943962
ISO Reg No. QACS-A-0623
US FDA Reg No. 3010077069

INDIA

Corporate Off. & Factory:
Kinra Park P.O.
Thrissur Dist., Kerala
India - 680 309
+91 480 2734989

JJ design & manufacture complete Surgical Instruments for - General Orthopaedic, Spine, Hand & Micro Surgery, Gynec & Obstetrics, Uro, E.N.T, Neuro, Thoracic & Cardiovascular, Plastic, General surgery.

Features and Benefits

- Genuine 'JJ' Brand
- Custom made instruments
- Made out of German Stainless Steel- SS410& SS 420
- Follows International Passivation standard - ASTM 967
- International hardness of product - 45 HRC & 50 HRC
- Universal Specifications & Quality
- 5 years replacement warranty against manufacturing defect
- Absolute Competitive price
- JJ welcomes innovative ideas and designs
- Gift packs on attractive terms

JJ's repairing & refurbishing program:

- Repairing & refurbishing of surgical instruments of any brand
- Reduce purchase budget drastically
- Annual Maintenance Contract (AMC)
- Round the clock service center
- Replace damaged instrument with new until rectification

JJ's Mobile Trade Expo

- Introducing mobile sales and service all over India
- Specially designed instruments on display
- Easy to select, exchange and experience
- Get technical advice from the experts
- Contact us for mobile trade expo schedules

www.myjjonline.com

sales@myjjonline.com, myjjonline@yahoo.com

UAE

Sales Office:
Suite E-74G-35
Hamriyah Free Zone
Sharjah, UAE
+971 50 7571388

OMAN

Service & Sales Office:
Suite 122, Bldg 579
Masabala Industrial Area
Muscat, Oman
+968 9905 3791

USA

Sales Office:
1000 Centre Greenway
Suite 208, Cary
North Carolina 27513
+1 919 467 5278



INDEX



	Page No.
SOFT TISSUE INSTRUMENTS	
Easy-Cut Scissors	1 & 4
Dissecting Forceps	5 & 6
Needle Holders	7 & 8
Skin Hooks	9
Hand Held Retractors	10
Self Retaining Retractors	11
Tuppers Hand Retractor	12
Artery Forceps	13
Scalpel Handles & Marking Pens	14



BONE INSTRUMENTS	
Bone Nibblers	15 & 16
Bone Cutters	17
Mallets	18
Osteotomes & Gouges	19 & 20
Bone Elevators	21 & 22
Bone Curettes	23
Dissectors	24
Bone Holding Forceps	25
Bone Spreaders	26
Wiring Cutters, Pullers & Pliers	27 & 28




TENDON INSTRUMENTS	
Tendon Interlacing Forceps & Retrievers	29 & 30
Tendon Strippers	31
Meyer Tendon Trimming Set	32
Recommended Hand Instrument Set	33 & 34



SMALL JOINT ARTHROSCOPY INSTRUMENTS	
Arthroscopic Punch Forceps	35
Arthroscopic Hooks	36
Arthroscopic Telescopes	37 & 38



MICRO INSTRUMENTS	
Micro Vessel Dilators	39 & 40
Micro Forceps	41 & 42
Micro Scissors	43 & 44
Micro Needle Holders	45
Acland Micro Vessel Clamps	46 & 47
Recommended Micro Instrument Sets	48
Instrument Containers	49
Sterile Packed Kirschner Wires	50
Digiband Sterile Digital Tourniquet	51



Iris Scissors

- 11.5cm/4½" straight
- Micro serrated cutting edge


HF7612-11



Iris Scissors

- 11.5cm/4½" curved
- Micro serrated cutting edge

HF7613-11



Strabismus Scissors

- 11.5cm/4½" straight
- Micro serrated cutting edge

HF7622-11



Strabismus Scissors

- 11.5cm/4½" curved
- Micro serrated cutting edge

HF7623-11



Stevens Tenotomy Scissors

- 11.5cm/4½" curved
- Micro serrated cutting edge


HF7633-11



Stevens-Stille Pattern Tenotomy Scissors

- 12.5cm/5" curved
- Micro serrated cutting edge
- Larger oversized rings

HF7633-12



Jameson-Werber Scissors

- 13.5cm/5¼" curved
- Micro serrated cutting edge

HF7643-13



Jameson Scissors

- 15cm/6" curved
- Micro serrated cutting edge

HF7643-15



Jameson Scissors

- 15cm/6" curved
- Strong pattern
- Micro serrated cutting edge

HF7645-15




SCISSORS

Slim Baby Metzenbaum Easy-Cut Scissors

- 14cm/5½" curved
- Micro serrated cutting edge


HF7661-14



Slim Metzenbaum Easy-Cut Scissors

- 18cm/7" curved
- Micro serrated cutting edge

HF7661-18



Baby Metzenbaum Easy-Cut Scissors

- 14cm/5½" curved
- Micro serrated cutting edge

HF7663-14



McIndoe / Metzenbaum Easy-Cut Scissors

- 18cm/7" curved
- Micro serrated cutting edge

HF7663-18



SCISSORS

Mayo Easy-Cut Scissors

- 14cm/5½" straight
- Micro serrated cutting edge


HF7652-14



Mayo Easy-Cut Scissors

- 14cm/5½" curved
- Micro serrated cutting edge

HF7653-14



Mayo Easy-Cut Scissors

- 17cm/6¾" straight
- Micro serrated cutting edge

HF7652-17



Mayo Easy-Cut Scissors

- 17cm/6¾" curved
- Micro serrated cutting edge

HF7653-17




DISSECTING FORCEPS

Gillies Dissecting Forceps

Toothed
15cm/6"

HF5160-15



McIndoes Dissecting Forceps

Non toothed
15cm/6"

HF5170-15



Debakey-Diethrich Tissue Forceps

1.5mm tips
15cm/6"

HF2015-15



Debakey Tissue Forceps

2mm tips
15cm/6"

HF2120-15



Bipolar Forceps

1mm tips
11cm/4½"

70-02-41



Bipolar Forceps

1mm tips
16cm/6¼" length

70-01-51



DISSECTING FORCEPS

Adson Forceps

Non toothed
12cm/4¾"


HF5025-12



Adson Forceps

Toothed
12cm/4¾"

HF5120-12



Micro Adsons Forceps

Non toothed
12cm/4¾"

HF5035-12



Micro Adson Forceps

Toothed
12cm/4¾"

HF5130-12



Adson Brown Forceps

7x7 teeth
12cm/4¾"


HF5140-12



Adson-DeBakey Forceps

2mm tips
12cm/4¾"

HF5180-12



NEEDLE HOLDERS

Micron Fine Needle Holder

12.5cm/5"

Diamond dusted jaws



HF1018-12



Derf Needle Holder

12cm/4¾"

Tungsten Carbide jaws



HF703DE12

Smooth



HF705DE12



Halsey Needle Holder

13cm/5"

Tungsten Carbide jaws



HF706.DE-13

Smooth



HF704.DE-13



Crile Wood Needle Holder

Tungsten Carbide jaws



HF712.DE-15

15cm/6"



HF734.DE-18

18cm/7"



HF734.DE-20

20cm/8"




NEEDLE HOLDERS

Nieverts Needle Holder (offset ring)

13cm/5¼"

Tungsten Carbide jaws



HF783.DE-13



Foster Gillies Combined Needle Holder/Scissors

13cm/5"

Tungsten Carbide jaws




HF775.SE-13



Gillies Combined Needle Holder/Scissors

16cm/6¼"

Tungsten Carbide jaws



HF775.DE-16



SKIN HOOKS

Kilner Detachable Skin Hook

15cm/6"

HF5075-01



Gillies Skin Hook

18cm/7"

HF5065-01

Small size



HA5065-02

Large size



Klienert-Kutz Skin Hook

16cm/6¼"

HF5055-01



Joseph Double Skin Hook

16cm/6¼"

HF5402-10



Blunt Nerve Hook

16cm/6¼"

HF5251-01



Blunt Crile Nerve Hook

14cm/5¾"

HF5200-14



Cushing Blunt Hook

19cm/7½"

3mm tip

HF9101-19



HAND HELD RETRACTORS

Catspaw Retractors

16cm/6¼"
Semi-sharp prongs

HF5252-03



Senn-Miller Retractor

16cm/6¼"
Double ended

HF5289-01 Sharp prongs

HF5289-02 Blunt prongs



Kilner-Miniature Pattern Catspaw Retractor

15cm/6"
Double ended

HF5286-01



Ragnell-Davies Retractor

14cm/5½"
Double ended

HF5287-14



Desmarres Saddle Retractor

16cm/6¼"
12mm width

HF5235-12



Mini-Langenbeck Retractor

16cm/6¼"

HF5242-10 10 x 6mm*

HF5242-17 17 x 5mm



*illustrated

Langenbeck Retractor

• 21.5cm/8½"



	Blade Size		
	depth	width	
SP5122-21	22mm	6mm	Small
SP5145-21	45mm	13mm	Medium
SP5165-21	65mm	25mm	Large

SELF RETAINING RETRACTORS

Alms Retractor

4 x 4 teeth blunt
7cm/2¾"
4mm depth prongs

HF5200-70



Baby Weitlaner S/R Retractor

2 x 3 blunt prongs
11.5cm/4½"

HF5701-10



Wests S/R Retractor

3 x 4 sharp prongs
14cm/5½"

HF4003-14



Weitlaner S/R Retractor

3 x 4 blunt prongs
14cm/5½"

HF4002-14



Finsen Small Wound Retractor

3 x 4 Teeth Blunt
7cm/2¾"

HF5705-10



Kaupp Wound Retractor

4 x 4 Teeth Blunt
10cm/4"



HF5710-14



HAND RETRACTOR

Tuppers Hand Retractor Set
Complete Set with container

HF-100/110




Tuppers Universal Hand Holder and Retractor Set enables the surgeon to fixate the hand in any position needed. The patient's hand is held absolutely firm and allows unhindered work under the microscope.

Ball chains of varying length with skin or tendon hooks fixate the hand and lower arm on the Hand Holder.

Due to the slots around the Hand Holder, retraction in every direction is possible.

Upward retraction is easily achieved with the elevated chain holders.

The special container is designed for easy handling and safe storage.



Set consisting of the following instruments:

Qty	Order Code	Description
1	TH-101	Hand Table
2	TH-102	Elevated chain holders
2	TH-103	Chain holder hooks
2	TH-104	Ball chains 32cm, with skin-hooks
4	TH-105	Ball chains 21cm, with skin-hooks
2	TH-106	Ball chains 26cm, for arm and wrist fixation
4	TH-107	Tendon hooks 5.0mm wide
4	TH-108	Sliding hooks for the rubber bands
1	TH-109	Set rubber bands 5 pieces

ARTERY FORCEPS

Micro Baby Mosquito Artery Forceps

Straight
10cm/4"

HF5742-10



Micro Baby Mosquito Artery Forceps

Curved
10cm/4"

HF5743-10



Halstead Mosquito Artery Forceps

Straight
12.5cm/5"

HF5752-12



Micro Mosquito Artery Forceps

Straight
12.5cm/5"

HF5742-12



Halstead Mosquito Artery Forceps

Curved
12.5cm/5"

HF5753-12

Micro Mosquito Artery Forceps

Curved
12.5cm/5"

HF5743-12

Barrons Scalpel Blade Handle No.3

13cm/5¼"

HF5001-13



Beaver Blade Holder

13cm/5¼"

HF5004-13



Eckoff Mapping Pen

14.5cm/5¼"

HF5002-14



Sommerlad Marking Pen

15cm/6"

HF5003-15



Stainless Steel Measuring Ruler

Code Length

HF5005-15

15cm/6"

HF5005-30

30cm/12"



Mitchels Trimmer

16cm/6¼"

HF5008-16



Castroviejo Measuring Caliper

7cm/2¾"

Measuring range 0 to 25mm

HF7334-10



Ink Pot

Stainless Steel

HA5006-99



BONE NIBBLERS

Friedman-Pearson Micro Bone Nibbler

Curved
15cm/6"

1mm Bite

HF5117-15



Kleinert-Kutz Bone Nibbler

Curved
14cm/4½"

2mm Bite

HF5118-50



Stellbrink Synovectomy Nibbler

Curved
17cm/6¾"

2mm Bite

HF5119-00



Beyer Bone Nibbler

Curved
17cm/6¾"

2mm Bite

HF5121-17



BONE NIBBLERS

Bohler Bone Nibbler
Curved double action
15cm/6"

3mm Bite


HF5120-50



Beyer Bone Nibbler
Curved double action
18cm/7"

4mm Bite


HF1002-18



Jansen Bone Nibbler
Curved double action
18cm/7"

5mm Bite

HF1003-18



BONE CUTTERS

Bohler Bone Cutters

Double action
15cm/6"

Straight

HF5160-50



Angled

HF5161-50



McIndoe Bone Cutters

Double action
18cm/7"

HF1028-18



Ruskin-Listons Bone Cutters

Double action
19cm/7½"

Straight

HF5170-19



Angled

HF5171-19



MALLETS

Toffee Mallet

19cm/7½"
Head size 22mm Ø
Weight 165g/6oz

HF5808-21



Mead Mallet

Interchangeable delrin discs
18cm/7"
Head size 26mm Ø
Weight 320g/11¼oz

HF5805-18



Ultra Fine Mallet

17cm/6¾"
Head size 20mm Ø
Weight 145g/5oz

HF5801-17



Mallet with Tufnol Handle

23cm/9"
Head size 35mm Ø
Weight 550g/20oz

HF5804-23



London Splint Mallet

18cm/7"
Head size 30mm Ø
Weight 450g/16.2oz

HF5802-18



Mini-Lambotte Osteotomes

12.5cm/5"

- HF5270-28 Set of 5 sizes
2mm to 8mm
- HF5270-02 2mm
- HF5270-03 3mm
- HF5270-04 4mm
- HF5270-06 6mm
- HF5270-08 8mm



Sheehan Osteotomes

15cm/6"

- HF5280-28 Set of 5 sizes
2mm to 10mm
- HF5280-02 2mm
- HF5280-04 4mm
- HF5280-06 6mm
- HF5280-08 8mm
- HF5280-10 10mm



Carroll Osteotomes

14.5cm/5¾"

- HF5290-03 3mm
- HF5290-06 6mm
- HF5290-10 10mm
- HF5290-13 13mm



Partsch Chisel

13cm/5½"

- HF5260-02 2mm
- HF5260-03 3mm
- HF5260-04 4mm
- HF5260-05 5mm
- HF5260-06 6mm
- HF5260-07 7mm
- HF5260-08 8mm



Ultra Fine Osteotomes (Titanium Nitrate Coated)

16cm/6¼"

HF5230-48 Set of 6 sizes
2mm to 12mm



HF5230-02 2mm

HF5230-04 4mm

HF5230-06 6mm

HF5230-08 8mm

HF5230-10 10mm

HF5230-12 12mm

Mini-Lexer Osteotomes with Tufnol Handle

18cm/7"

HF5250-48 Set of 5 sizes
2mm to 10mm



HF5250-02 2mm

HF5250-04 4mm

HF5250-06 6mm

HF5250-08 8mm

HF5250-10 10mm

Mini-Lexer Gouge with Tufnol Handle

18cm/7"

HF5240-04 4mm

HF5240-06 6mm

HF5240-08 8mm

HF5240-10 10mm



Bone Tampers

16cm/6½"

HF5004-02 2mm #

HF5004-03 3mm #

HA5004-04 4mm 

HA5004-05 5mm 

HA5004-06 6mm 

HA5004-08 8mm 



BONE ELEVATORS

Adson Periosteal Elevator

16cm/6¼"
Curved blunt



Width

HF3055-16	5mm
HF3007-16	7mm

Pennybacker Periosteal Elevator

16cm/6¼"
Curved sharp



Width

HF3063-16	4mm
HF3066-16	6mm
HF3068-16	8mm

Bone Raspatory with Tufnol Handle

19cm/7½"



Straight

Width


HF7016-19	6mm
-----------	-----



Curved

Width

HF7023-19	3mm
HF7026-19	6mm



Farabeuf Rugine Elevator

15cm/6"



Straight

Width

HF3001-15	12mm
-----------	------

Curved*

Width

HF3002-15	12mm
-----------	------

*illustrated

BONE ELEVATORS

Mini Bone Lever

14cm/5½"
2-8mm width

HF5830-14



Hohmann Style Mini-Bone Lever

16cm/6¼"
6mm width

HF5831-16



Phalangeal Lever

15cm/6"
5mm width

HF5832-16



Buck-Gramcko Lever

15cm/6"
7.5mm width

HF5833-15



Mini-Cobra Bone Lever

16cm/6¼"
15mm width

HF5834-16



McGlamry Metatarsal Elevator

17cm/6½"
Sharp



HF5835-10 9mm Width

HF5835-12 11mm Width

HF5835-14 13mm Width

HF5835-16 15mm Width

HF5835-19 19mm Width

BONE CURETTES & RASPS

Martini Mini Curette

13.5cm/5¼"

HF5412-13



Jansen Curette

16cm/6¼"

HF5432-15



Bruns Curette

17cm/6¾"



Cup size 0000

HF5414-17



Cup size 000

HF5413-17



Cup size 00

HF5412-17



Cup size 0

HF5411-17



Curette with Tufnol Handle

19cm/7½"



Cup size 0000

HF5404-19



Cup size 000

HF5403-19



Cup size 00

HF5402-19



Cup size 0

HF5401-19



Bone Rasp with Finger Ring

15cm/6"

HF5382-01



DISSECTORS

Swedish Dissector

18cm/7"



HF9201-18

Howarth Elevator

21.5cm/8½"



HF9207-21

MacDonald Elevator/Dissector



Size

HF7105-191

8cm/7"

HF9206-23

23cm/9"

Watson Cheyne Dissector



Size

HF7203-18

1 8cm/7"

HF9203-23

23cm/9"

HF9203-25

25.5cm/10"

Freers Elevator/Dissector

18.5cm/7¼"



HF7602-18

BONE HOLDING FORCEPS

Bone Holding Forceps

14cm/5½"

HF5021-14



Phalanx Holding Forceps

14cm/5½"


HF5031-14



Langenbeck Bone Holding Forceps

20cm/8"


HF5120-05



Bone Reposition Forceps for small fragments

13.5cm/6¼"

HF5041-13




Bone Repositioning Forceps with K-wire guide

14cm/5½"

HF5051-09 for 0.9mm wire

HF5051-12 for 1.2mm wire

HF5051-16 for 1.6mm wire



These Bone Repositioning Forceps include a wire guide allowing the K-wire to be accurately implanted in the Bone. We offer sizes to take small, medium and large Kirschner wires.

BONE SPREADERS & CAMPS

Scaphoid Spreader - Fine Pattern

14cm/5½"

1 x 2 prongs

6mm

HF5131-14



Scaphoid Spreader - Standard Pattern

14cm/5½"

6mm

HF5134-14



Baby Inge Lamina Spreader

18cm/7"

7mm
Toothed



HF5134-18

7mm
Smooth



HF5134-18S

7mm
Cross Serrated



HF5134-18X



Ikuta Bone Clamp

12cm/4¾"

Toothed 2x2 jaws

With ratchet




HF5055-12

Wide jaw*
8mm width



HF5050-12

Narrow jaw
4mm width



*illustrated

Muller Claus Wire Cutters

- 16cm/6"
- Tungsten carbide edges
- Cuts wire up to 1.6mm Ø

8-810E16



Muller Claus Wire Cutters

- Tungsten carbide edges

8-810E18

- 18cm/7" length (*illustrated)
- Cuts wire up to 2mm Ø

8-800E22

- 22cm/8¾" length
- Cuts wire up to 3mm Ø



End Wire Cutters


- 18cm/7" length
- Tungsten carbide edges
- Cuts wire up to 1.7mm Ø

8-816E18




Muller Claus Wire Pullers

- 18cm/7" length




2mm Jaws

8-818E18



4mm Jaws

8-820E18



6mm Jaws

8-821E18



Slim Nose Wire Twisters

- 19cm/7½"

HF7645-19



Needle Nosed Locking Pliers



Length

200-514-220	22cm/8¾"	Small
200-515-250	25cm/9¾"	Medium
200-516-300	30cm/11¾"	Large

- A very useful instrument for orthopaedic surgery allowing the surgeon to securely clamp on to wires and implants.

Combined Wire Cutters & Pliers

- 18cm/7"
- Suitable for wires up to 1.6mm

200-512-185



Orthopaedic Locking Grips

- Hardened jaws

Length

200-511-180	18cm/7"	Small
200-511-200	20cm/8"	Medium
200-511-240	24cm/9½"	Large



Tendon Interlacing Forceps

- 15cm/6"




Straight

HF5254-15



Curved

HF5255-15



The pointed tip separates the Tendon without causing splits. The jaws grasp and secure a Tendon Stump pulling it back through the tendon body.

Carroll Tendon Retriever

- Curved
- 11cm/4¼"


HF5679-11




Kleinert-Kutz Tendon Retriever

- 21cm/8¼"
- 2mm jaws
- Flexible shaft


HF5814-02



1 x 2 Teeth




Anderson-Brand Tendon Retriever



15cm/6"		20cm/8"	
HF5810-15	Straight	HF5810-20	Straight
HF5811-15	Curved	HF5811-20	Curved

Buck-Gramcko Tendon Retriever

- 20cm/8"
- 2 x 3 teeth on serrated jaws




HF5824-20	Straight	HF5825-20	Curved
-----------	----------	-----------	--------



Steinmann Tendon Holding Forceps

- 13cm/5 1/8"
- with adjustable grub screw

HF5864-13



Teeth engage in opposing jaws




TENDON SURGERY


Tendon Strippers


• 28cm/11"



Closed Ring Style

HF5811-03 3mm inside diam. 


HF5811-04 4mm inside diam. 

HF5811-05 5mm inside diam. 


HF5811-06 6mm inside diam. 

Slotted Ring Style

• No requirement to cut the Tendon

HF5809-03 3mm inside diam. 

HF5809-04 4mm inside diam. 

HF5809-05 5mm inside diam. 

HF5809-06 6mm inside diam. 

Fischer Tendon Stripper

• 31cm/12¼"



Diameter 4mm



HF5812-04



Larson Tendon Holding Forceps

• 10cm/4"



HF5810-02 2mm 

HF5810-03 3mm 

HF5810-04 4mm 

HF5810-05 5mm 

TENDON SURGERY

Meyer Nerve/Tendon Holding & Trimming Instrument Set

NHS-1

Set of Six Forceps with Container and Blade Holders.
CBS - 35 Box of Blades to fit Meyer Set.

Set consisting of the following instruments:

Qty	Order Code	Description
1	SHF-2.0	Holding Forceps
1	SHF-3.0	Holding Forceps
1	SHF-4.0	Holding Forceps
1	SHF-5.0	Holding Forceps
1	SHF-6.0	Holding Forceps
1	SHF-7.0	Holding Forceps
1	SHS-12	Bladeholder
1	HFN-913	Sterilization Box



BHS-12
Blade Holder



Sterile Blades

- Box of 12 straight blades

HFS-35

This set was developed to allow the surgeon to hold securely and trim the nerve ends with the correct size holding forceps avoiding trauma to the nerve. There is minimal loss of nerve end with this technique.

Holding Forceps



SHF-0.5 Ø 0.5mm



SHF-1.0 Ø 1.0mm



SHF-1.5 Ø 1.5mm



SHF-2.0 Ø 2.0mm



SHF-2.5 Ø 2.5mm



SHF-3.0 Ø 3.0mm

SHF-3.5 Ø 3.5mm

SHF-4.0 Ø 4.0mm

SHF-5.0 Ø 5.0mm

SHF-6.0 Ø 6.0mm

SHF-7.0 Ø 7.0mm

HAND SURGERY INSTRUMENTS SETS

The JJ Hand Set gives the consultant a range of instrumentation for both Soft Tissue and Bone Surgery, presented in an instrument container.

The set can be ordered as a Basic Set or the extended set with further Bone and Tendon Instruments.





JJ Basic Hand Set

HF5000-00

Set consisting of the following instruments:

Qty	Order Code	Description
1	E4026/5D	Sterilising Container with double rack
1	EA7633-11	Stevens Tenotomy Scissors Curved 4½"
1	703 DE-12	Derf TC Needle Holder 4¾"
1	HF5118-50	Kleinert Kutz Bone Nibbler 5½"
1	HF5200-70	Alm S/R Retractor 3"
1	HF5701-10	Weitlaner S/R Retractor 2 x 3 Blunt 4¼"
1	HF5130-12	Micro Adson Forceps 1 x 2 Teeth
1	HF5035-12	Micro Adson Forceps Non Toothed
1	HF5382-01	Bone Rasp Straight with Finger Ring
2	HF5055-01	Kleinert Kutz Skin Hooks
2	HF5402-10	Joseph Double Hooks 10mm
2	HF5251-01	Single Blunt Hook
2	HF5235-12	Saddle Hooks 12mm
2	HF5242-10	Mini Langenbeck Retractors
2	HF5252-03	Catspaw Retractors 3 - Prong
5	HF5270-28	Mini-Lambotte Osteotomes 2, 3, 4, 6 & 8mm
1	HF5432-10	Jansen Mini Bone Curette
2	HF5830-14	Mini Bone Levers
1	HF7602-18	Freer Dissector
1	HF5808-21	Toffee Hammer



JJ Extended Hand Set

HA5000-05


Extended Hand Set includes the above plus the following instruments:

Qty	Order Code	Description
1	EA7643-13	Jameson-Werber Scissors Curved 5¼"
1	70-02-41	Bipolar Forceps 4"
1	37-01-35L	Bipolar Lead - Banana Plug Fitting
1	HF5120-50	Bohler Bone Nibbler Curved
1	HF5161-50	Bohler Bone Cutter Angled
1	8-810 E-18	Muller-Claus Wire Cutters
1	8-818 E-18	Wire Twisters Muller-Claus
1	HF5679-11	Carroll Tendon Forceps
1	HF5254-15	Tendon Braiding Forceps Straight
1	HF5255-15	Tendon Braiding Forceps Curved

Vessel Dilator Angled

- 11cm/4¼" length
- Tips 0.2mm Ø

HFD-5a.2



Vessel Dilator Angled

- 11cm/4¼" length
- Tips 0.3mm Ø

HFD-5a.3



Vessel Dilator

- 13.5cm/5³/₈" length
- Tips 0.2mm Ø

HFS-3d.2



Vessel Dilator

- 13.5cm/5³/₈" length
- Curved tips 0.2mm Ø

HFCL-3d.2



Vessel Dilator Long

- 18cm/7" length
- Tips 0.3mm Ø


HF-3-18d.3



Vessel Dilator Long

- 18cm/7" length
- Curved tips 0.2mm Ø

HFCL-3-18d.2



JJ Vessel Dilators have Mirror polished parallel tips for interluminal vessel dilation. The entire area is opened to the same diameter and the pressure is distributed evenly. This ensures that unnecessary damage to the vessel is avoided.

Round Bodied Vessel Dilator

- 15cm/6" length
- Straight tips 0.3mm Ø
- Balanced forceps

HFS-15 RM-8d.3



Round Bodied Vessel Dilator

- 15cm/6" length
- Curved tips 0.2mm Ø
- Balanced forceps

HFC-15 RM-8d.2



Round Bodied Vessel Dilator

- 18cm/7" length
- Straight tips 0.3mm Ø
- Balanced forceps

HFS-18 RM-8d.3



Round Bodied Vessel Dilator

- 18cm/7" length
- Curved tips 0.3mm Ø
- Balanced forceps

HFC-18 RM-8d.3



Jewellers Forceps

- 11cm/4¼" length
- Tips 0.3mm

HF-5



Jewellers Forceps Long

- 15cm/6" length
- Tips 0.3mm

HFL-5



Jewellers Forceps

- 12cm/4¾" length
- Tips 0.3mm

HF-3



Jewellers Forceps

- 13.5cm/5¾" length
- Tips 0.3mm

HFL-3



Pierse Forceps

- 12cm/4¾" length
- Tips 0.3mm

HFP-3



Pierse Forceps have a micro hole in the tips which securely grip tissue with minimal trauma



Micro Adson Fixation Forceps

- 12cm/4¾" length
- 1 x 2 micro teeth
- Tips 0.3mm

HF-1



Round Bodied Forceps

- 15cm/6" length
- Straight, with platform
- Tips 0.3mm Ø
- Balanced forceps

HFS-15 RM8



Round Bodied Forceps

- 15cm/6" length
- Curved, with platform
- Tips 0.3mm Ø
- Balanced forceps

HFC-15 RM8



Round Bodied Angled Forceps

- 15cm/6" length
- Angled tips 0.2mm Ø
- Balanced forceps

HFAS-15 RM-8



Round Bodied Forceps

- 18cm/7" length
- Straight, with platform
- Tips 0.5mm Ø
- Balanced forceps

HFS-18 RM8



Round Bodied Forceps

- 18cm/7" length
- Curved, with platform
- Tips 0.5mm Ø
- Balanced forceps

HFC-18 RM8



Round Bodied Angled Forceps

- 18cm/7" length
- Angled tips 0.2mm Ø
- Balanced forceps

HFS-18 RM8




Straight Scissors Flat Handle

- 15cm/6" length
- Blades 19mm long

HFS-15



Adventitia Scissors with pointed tips



SAS-15T

Available with Micro Serrations for nerve cutting -Straight


Curved Scissors Flat Handle

- 15cm/6" length
- Blades 19mm long

HFC-15



Dissecting Scissors with rounded tips



SDC-15T

Available with Micro Serrations for nerve cutting -Curved


Straight Scissors Flat Handle

- 15cm/6" length
- Blades 15mm long

HFS-15S



Adventitia Scissors with pointed tips



Curved Scissors Flat Handle

- 15cm/6" length
- Blades 15mm long

HFC-15S



Dissecting Scissors with rounded tips



Straight Scissors Round Handle

15cm/6" SAS-15 R8


- Blades 12mm long

18cm/7" SAS-18 R8

- Blades 14mm long



Adventitia scissors with pointed tips



Curved Scissors Round Handle

15cm/6" HFC-15 R8

- Blades 12mm long

18cm/7" HFC-18 R8

- Blades 14mm long



Dissecting scissors with rounded tips



Blondeel Scissors Round Handle

- Blades 15mm long
- Curved blades

15cm/6" HFC-15R-8-2R

18cm/7" HFC-18R-8-2R



These Scissors were designed by professor Phillip Blondeel for use in breast reconstruction utilising a perforator flap. A combination of ring handle design for preparation of fine structures, the round handle micro instrument handle allows the user to perform fine micro surgery.

Barraquer Micro Needle Holder

- Round bodied
- 0.4mm Jaws



Without Lock	Straight	Curved
13cm/5"	BS-13-8	B-13-8
15cm/6"	BS-15-8	B-15-8
18cm/7"	BS-18-8	B-18-8

Castroviejo Micro Needle Holder

- Flat handled
- 0.4mm jaws




Without Lock	Curved
14cm/5 1/2"	C-14
18cm/7"	C-18

Double Action Micro Needle Holder

- Round bodied
- 0.4mm jaws



Without Lock	Curved
15cm/6"	NH-15 RM-8
18cm/7"	NH-18 RM-8




"A" or "V" Pattern?

Throughout the series, the "V" pattern is the all-purpose clamp, commonly used on all veins, and almost all arteries. The "A" pattern, available in all sizes, is designed for exceptionally thick-walled vessels and slippery arteries. The slight incurve of the tips of the "A" pattern jaw provides extra security against slippage.

- Light and compact design
- Extremely durable
- A unique gripping surface
- Smooth sliding bar action on all approximator clamps
- Clamp profile with parallel sided jaws
- Matt or black finish

















Acland
Clamp Container
HF4400-00




Double
Clamps

Single
Clamps

Ordering details

SINGLE CLAMPS	DOUBLE APPROXIMATING CLAMPS	DOUBLE CLAMPS WITH FRAME
00396 V 00397 A  CLAMP SIZE B - 1 0.4 - 1.0mm	00414 V 00415 A  CLAMP SIZE ABB - 11 0.4 - 1.0mm	00408 V 00409 A  CLAMP SIZE ABB - 1 0.4 - 1.0mm
00398 V 00399 A  B - 2 0.6 - 1.5mm	00416 V 00417 A  ABB - 22 0.6 - 1.5mm	00410 V 00411 A  ABB - 2 0.6 - 1.5mm
00400 V 00401 A  B - 3 1.0 - 2.25mm	00418 V 00419 A  ABB - 33 1.0 - 2.25mm	00412 V 00413 A  ABB - 3 1.0 - 2.25mm
00325  HD - S 1.5 - 3.5mm	00326  HD - D 1.5 - 3.5mm	 <p>Clamp Box CB-1</p> <p>Clamp Applicator CAF-4</p>
00286  RD - S 2.0 - 5.0mm	00287  RD - D 2.0 - 5.0mm	



- Sterile, single use, disposable
- Atraumatic jaws
- Colour coded for better differentiation
- Consistent and reliable closing pressure
- Designed for use on veins and arteries
- Light weight
- Blue, black and purple finish provides excellent contrast in the field
- Radiopaque
- Eliminates cleaning and processing problems
- Economical

Single Clamps packed in boxes of 24



P - 1
0.4 - 1mm
Vessels

HF70005



P - 2
0.6 - 1.5mm
Vessels

HF70003



P - 3
1.0 - 2.5mm
Vessels

HF70007

Double Clamps packed in boxes of 12



APP - 11
0.4 - 1mm
Vessels

HF70006



APP - 22
0.6 - 1.5mm
Vessels

HF70004



APP - 33
1.0 - 2.5mm
Vessels

HF70008


VISIBILITY BACKGROUND MATERIAL

TO ASSIST IN VISUALISING MICRO VESSELS DURING SURGERY



Application: For use in Micro Surgical Anastomosis of Micro Vessels under a Microscope to allow the surgeon to better visualize the detail of the vessels, being closed.

The Visibility Material can be used either with the Grid side up to give the Surgeon a scale of the vessel alternatively the reverse Plain side can also be used.

The material has a none reflective surface eliminating distractions from theatre lighting and has a non-tacky finish to aid mobilization of the Vessel.



1mm Grid one side & matt finish on reverse. Size: 50mm x 25mm



Ordering details:

COLOUR	ordering code
Yellow	VB2
Green	VB3
Blue	VB4

- X-Ray opaque silicon material
- 1mm Grid for accurate measurement of vessels
- Latex Free material
- Sterile Packed in boxes of 12 pieces

JJ offer both a Basic Micro Surgery Instrument Set, equipping the surgeon with the basic Instruments for most micro surgery procedures and a reduced set suitable for Nerve and Tendon Repair.

JJ Basic Micro Instrument Set

BMS-1

Set consisting of the following instruments:



Qty	Order Code	Description
1	E.4026/4	Micro Instrument Container
1	SAS-15R8	Micro Scissors R/B 15cm Straight
1	SDC-15R8	Micro Scissors R/B 15cm Curved
1	BS-15-8	Micro Needle Holder R/B 15cm Straight
1	B-15-8	Micro Needle Holder R/B 15cm Curved
1	FRS-15 RM8	Suture Tying Forceps R/B 15cm Straight
1	FRC-15 RM8	Suture Tying Forceps R/B 15cm Curved
2	JF-5	Jewellers Forceps No.5 11cm
1	D5a	Micro Vessel Dilator 11cm
1	CAF-4	Clamp Applying Forceps - No Lock
1	CB-1	Clamp Box
1pr	00398V	Single Acland Clamps B-2V
1pr	00400V	Single Acland Clamps B-3V
1pr	00325V	Single Acland Clamps HD-S

Micro Instrument Set for Nerve Repair Surgery

MNS-1

Set consisting of the following instruments:



Qty	Order Code	Description
1	E.4026/1	Micro Instrument Container
1	B-15-8	Micro Needle Holder R/B 15cm Curved
1	SAS-15T	Micro Scissors Flat Handled 15cm Straight Serrated Blade
2	JF-5	Jewellers Forceps No.5 11cm




SINGLE TRAY CONTAINERS

Internal Dimensions			Code
A	B	C	
200 x	140 x	50mm	E 4026/1
280 x	195 x	50mm	E 4026/4
265 x	255 x	50mm	E 4026/5
380 x	260 x	50mm	E 4026/7

E 4026/5

Standard Container illustrated

One of our popular cases for Micro Instruments



DOUBLE TRAY CONTAINERS

Internal Dimensions			Code
A	B	C	
280 x	195 x	65mm	E 4026/4 D
265 x	255 x	65mm	E 4026/5 D
380 x	260 x	65mm	E 4026/7 D

E 4026/5D

Double Tray Container illustrated

Two removable perforated trays



Acland Clamp Container

HF4400-00

With individual compartments for various clamp sizes



International Instruments

Designers & Manufacturers of Surgical Instruments

Corporate Office & Factory

337/N, Amala Antonian Plaza
Kinfra Park P.O., Thrissur (Dt.)
Kerala State, India - 680 309
Tel : +91 480 273 5959 / 5858

USA Office

12916 Edsel Dr.
Raleigh NC 27613
USA
Tel : +1 919 467 5278

Middle East Office

211 Bait Al Ali Building
Al-Khuwair, Muscat
Oman
Tel : +968 99053791

sales@myjjonline.com, myjjonline@yahoo.com

www.myjjonline.com