




International Instruments

Designers & Manufacturers of Surgical Instruments

Laparoscopy Instruments



Orthopaedic | Spine | Hand & Micro
Plastic | Maxillofacial | Laparoscopy
General | Gynec & Obstetrics | Uro
E.N.T | Cardiovascular | Thoracic | Neuro

US FDA, CE & ISO REGD.

Kinfra Park P.O., Koratty, Thrissur District, Kerala, India 680309 | +91 480 2735959, 2735858, +91 94970 77555
sales@myjjonline.com | myjjonline@yahoo.com

www.myjjonline.com

INDIA • USA • OMAN • UAE


Trocars

LCT105



Multi Function 10.5 mm

LCT55



Multi Function 5.5 mm

LCT35



Multi Function 3.5 mm

LC011



Flip Top Plastic

LCTN105



New Multi Function 10.5 mm

LCTN55




New Multi Function 5.5 mm

LCF15



Flap 15 mm

LCF105




Flap 10.5 mm

LCF07



Flap 7 mm

LCF55




Flap 5.5 mm

Safety Trocars & Veress Needles


LCST10




LCHT10



LCST07



LCHT05



LCST05



LCVN120



LCVN140



LCVN100



LCVN160



LCDT05



Suctions & Protection Tube

LCTR10



Suction Trumpet 10

LCTH10



Suction Thumb Control 10

LCTR05



Suction Trumpet 5 mm

LCTH05



Suction Thumb Control 5 mm

LCPR10



Scope Protection Tube 10 mm

LCPR05



Scope Protection Tube 5 mm

Reducers & Insert Sheath

LCR20-10



20 / 10

LCR15-10



15 / 10

LCR10-07



10 / 7

LCR07-05



7 / 5

LCR05-03



5 / 3

LCIS01

Insert Sheath



Grasper Handles

LCGH1



Handle Plain (1)

LCGH2



Handle with Ratchet (2)

LCGH3



Handle Plain (3)

LCGH4



Handle with Ratchet (4)

Grasper Handles

LCGH5



Handle Plain (5)

LCGH6



Handle with Ratchet (6)

LCGH7



Handle with Ratchet (7)

Scissors & Graspers

LCSG1



D/A Curved Scissor

LCSG2



Maryland Flat Long

LCSG3



Bowl Long

LCSG4



D/A Straight Scissor

LCSG5



Maryland Long

LCSG6



Bowl Medium

LCSG7



D/A Short Scissor

LCSG8



Maryland Medium

LCSG9



Bowl Short

LCSG10



Micro Scissor

LCSG11



Maryland Short

LCSG12



Serrated

LCSG13



S/A Scissor

LCSG14



Babcock Long

LCSG15



Universal

LCSG16



Hook Scissor

LCSG17



Babcock Short

LCSG18



Serrated Short

LCSG19



Claw 2 x 3

LCSG20



Tinaculum

LCSG21



Oval Bowl

LCSG22



Claw 2 x 4

LCSG23



Debeky

LCSG24



Tube Holding

Scissors & Graspers

LCSG25



Allis

LCSG28



Allis Toothed

LCSG30



Mixer 90 deg

LCSG33



Mixer 45 deg

LCSG36



Biopsy

LCSG26



Endo Clinch Long

LCSG31



Endo Clinch Short

LCSG34



1 x 2 Rat Tooth

LCSG27



Alligator

LCSG29



Str Grasping

LCSG32



Maxi

LCSG35



Fundus

LCSG37



Dolphin

3 mm Inserts

LCSG38



Mixer 90 deg

LCSG39



Curved Scissor

LCSG40



Alligator

LCSG41



Maryland

LCSG42



Oval Bowl

LCSG43



Bowl

LCSG44



1 x 2 Tooth

LCSG45




Babcock

Needle Holder, Nande & TOT Needle

LCNH1

Storz Needle Holder



LCNH2

Ethicon Needle Holder



LCNH5



Straight Jaw

LCNH3



Self Riding Jaw

LCNH6



Parrot Jaw

LCNH4



Curved Jaw

LCNH7



1 / 2 Tooth Jaw

LCNNR

Nande Needle Right



Nande Needle Left

LCNNL




Nande Needle Curved

LCNNC



LCPCN

Port Closure Needle



LCTNR



TOT Needle R & L

LCTNL



Clip Applicators

LCCA1

Clip Applicator S/A - Adjustable Jaws



LCHCA5

Hemo Lock Clip Applicator 5 mm



LCCA5

Clip Applicator 5 mm




LCCA3

Clip Applicator 3 in 1 Jaws



LCHCA

Hemo Lock Clip Applicator



LCCA9


Clip Applicator 90 deg



Fan Retractors & Gold Finger

LCLR5

Liver Retractor 5 mm



LCLR10

Liver Retractor 10 mm




LCGF5

Gold Finger 90 deg 5 mm




LCFR10

Fan Retractor Elevated 10 mm, 5 Blades




LCFR5

Fan Retractor Elevated 5 mm, 5 Blades



LCFR53

Fan Retractor Elevated 5 mm, 3 Blades



LCFR9

Fan Retractor Elevated 90 deg 10 mm, 2 Blades



Bi-Polars

LCBP1 Bi - Polar Double Action

LCBP11 D/A Maryland

LCBP12 D/A Ruby

LCBP13 D/A Bowl

LCBP2 Bi - Polar Single Action Detachable

LCBP21 S/A Maryland

LCBP22 S/A Ruby

LCBP23 S/A Bowl

LCBPX Bi - Polar 'X' Handle

LCBPXS Spiral

LCBPXF Flat

LCBPU Bi - Polar 'U' Handle


LCBPUS Spiral

LCBPUF Flat

Bi-Clamp & Sealers


LCBCMN

Bi - Clamp Metal New




LCBCM

Bi - Clamp Metal




LCBCM2

Bi - Clamp Metal (2)




LCVS5

Vessel Sealer 5 mm



LCBFD

Bipolar Forceps Detachable



Uterine Manipulator


LCUMM

Uterine Manipulator Mangeshikar




LCUMD

Uterine Manipulator with Die Injector



LCUMS

Uterine Manipulator Storz



LCUMV

Uterine Manipulator V Type



Uterine Manipulator & Vaginal Tube

LCUSM



Sukadhya Manipulator

LCUMC



Uterine Manipulator Cup Set of 3

LCUMA



Marwa Manipulator

LCVT



Vaginal Tube

Claw Forceps, Mini Lap & Die Injectors

LCCF10

Claw Forcep 10 mm



LCS10

Spoon 10 mm



LCML

Mini Lap



LCMLT

Traumatic



LCMLNT

Non Traumatic



LCDI1

Die Injector 1



LCDI2


Die Injector 2



Working Element and & Urology

LCWEDS

Working Element Double Stem



LCDSDC

Diagonistic Sheath Double Channel



LCDSSC

Diagonistic Sheath Single Channel



LCTIOS

Turp Inner & Outer Sheath



LCOS

Operative Sheath



LCOSB

Operative Sheath Batuchi Type
2.9 & 4.00 mm Scope




Urology


LCWEA



LCWEP



LCSC



Cables

LCC





International Instruments

Designers & Manufacturers of Surgical Instruments

US FDA, CE & ISO REGD.

ABOUT US

JJ International Instruments - A leader in designing and manufacturing of high quality stainless steel surgical instruments in India. With a vast experience of 46 years in the industry, JJ provides unmatched service and unbeaten quality to the worldwide market. Corporate office and manufacturing unit with state of art facilities is located in south Indian state, Kerala - The renowned 'God's Own Country'.

CONTINUES SUCCESS...!

JJ is now worldwide since exhibited in AAOS 2013 with tremendous response from both domestic and international surgeons. JJ now participates several tradeshows in USA and more than 30 tradeshows all over India yearly.

" International Quality @ Indian Price "

CE Reg No. 2P191227.IIWT46
ISO 13485 Reg No. EUAC/MD/1014
ISO 9001 Reg No. EUAC/QMS/1065
US FDA Reg No. 3010077069

INDIA

Corporate Off. & Factory:
Kinfra Park P.O.
Thrissur Dist., Kerala
India - 680 309
+91 480 2735959

UAE

Sales Office:
Suite E-74G-35
Hamriyah Free Zone
Sharjah, UAE
+971 50 7571358

OMAN

Service & Sales Office:
Suite 122, Bldg 579
Maabela Industrial Area
Muscat, Oman
+968 9905 3791

USA

Sales Office:
1000 Centre Greenway
Suite 200, Cary
North Carolina 27513
+1 919 467 5278

JJ design & manufacture complete Surgical Instruments for - General Orthopaedic, Spine, Hand & Micro Surgery, Gynec & Obstetrics, Uro, E.N.T, Neuro, Thoracic & Cardiovascular, Plastic, General surgery.

Features and Benefits

- Genuine 'JJ' Brand
- Custom made instruments
- Made out of German Stainless Steel- SS 410 & SS 420
- Follows International Passivation standard - ASTM 967
- International hardness of product - 45 HRC & 50 HRC
- Universal Specifications & Quality
- 5 years replacement warranty against manufacturing defect
- Absolute Competitive price
- JJ welcomes innovative ideas and designs
- Gift packs on attractive terms

JJ's repairing & refurbishing program:

- Repairing & refurbishing of surgical instruments of any brand
- Reduce purchase budget drastically
- Annual Maintenance Contract (AMC)
- Round the clock service center
- Replace damaged instrument with new until rectification

JJ's Mobile Trade Expo

- Introducing mobile sales and service all over India
- Specially designed instruments on display
- Easy to select, exchange and experience
- Get technical advice from the experts
- Contact us for mobile trade expo schedules

www.myjjonline.com

sales@myjjonline.com, myjjonline@yahoo.com

