




International Instruments

Designers & Manufacturers of Surgical Instruments

USFDA, CE, ISO & CDSCO REGD.



www.myjjonline.com

sales@myjjonline.com, myjjonline@yahoo.com

INDIA

Corporate Off. & Factory:
Kinfra Park P.O.
Thrissur Dist., Kerala
India - 680 309
+91 480 2735959

UAE

Sales Office:
AL Smookh Building
Umm Al Quwain,
UAE
+971 50 7571358


OMAN

Service & Sales Office:
Suite 122, Bldg 579
Maabela Industrial Area
Muscat, Oman
+968 9905 3791

USA

Sales Office:
1000 Centre Greenway
Suite 200, Cary
North Carolina 27513
+1 919 467 5278

JJ Quality Certificates



The annual registration for the following medical device establishment has been successfully completed for 2024:

Registration Number: 3010077069
 Owner Operator Number: 10043535
INTERNATIONAL INSTRUMENTS
 Kinfra Park PO
 Koratty
 Thrissur, Kerala 680308
 INDIA


If you do not see a registration number assigned to the establishment and your establishment previously had one, please send an email to regist@cdrh.fda.gov and include the registration number you believe is assigned to your establishment. We will review and determine if a duplicate registration has been created for your establishment.

Your registration is valid until December 31, 2024. Registration for 2025 will be conducted between October 1 and December 31, 2024.

Please note that registering your device facility and listing your devices does not, in any way, constitute FDA approval of your facility or your devices.

Should you have any questions, please send an e-mail to the CDRH Registration and Listing Helpdesk at regist@cdrh.fda.gov.

CDRH Registration and Listing Helpdesk
 Imports & Registration and Listing Team
 Division 2 Establishment Support
 Office of Regulatory Programs
 Office of Product Evaluation and Quality
 Center for Devices and Radiological Health
 U.S. Food and Drug Administration
 Tel: 301-796-7400, Option 1
 Email: regist@cdrh.fda.gov




Certificate of Compliance

CE

We hereby declare that the technical file of product complied with the requirement of Medical Device Regulation (Council Regulation 2017/745).

Certificate No.: **CE-5862**

Manufacturer
 Name : **INTERNATIONAL INSTRUMENTS** 
 Address : **PLOT NO.XIX-350, KINFRA PARK POST-680309**
THRISSUR KERALA, INDIA

Products : **(AS PER APPENDIX OF D)**
 Complies with the requirements applicable to it

The quality system file has been assessed, approved and is subject to continuous surveillance according to the Medical Device Regulation (Council Regulation 2017/745)

This certificate is issued under the following conditions:

1. It applies only to the quality system maintained in the manufacture of above referenced models and it does not substitute the design or type examination procedures, if required.
2. The certificate remains valid until the manufacturing conditions or the quality systems are changed.
3. The certificate validity is conditioned by positive results or surveillance audits.

The CE mark as shown above can be used, under the responsibility of the manufacturer, after completion of an EC Declaration of conformity and compliance with all relevant EC Directives. The statement is based on a single evaluation of one sample of above mentioned product. It does not imply an assessment of the whole production.

Validity of this certificate can be verified at www.ukcertifications.org.uk/verify

Date of Certification	20th July 2024
1 st Surveillance Audit Due	19th July 2025
2 nd Surveillance Audit Due	19th July 2026
Certificate Expiry (subject to the company maintaining its system to the required standard)	19th July 2027


 Authorised Signatory



Page 1 of 2

This certificate is the property of UK Certification & Inspection Limited and shall be returned immediately on request.
 72-75 Shelton Street, Covent Garden, London, WC2E 8JF, United Kingdom
 Website: www.ukcertifications.org.uk, email: info@ukcertifications.org.uk
 Company No: 11847631

| evaluate | assess | certify |  assessors of quality

Certificate No. **MD-QMS/91/R/N/2643**

This is to certify that

INTERNATIONAL INSTRUMENTS

Kinfra Park P.O., Koratty, Thrissur - 680 309, Kerala, India

has been found to conform to the requirements of
 Medical Devices - Quality Management System Standard

ISO 13485:2016

This certificate is valid for the following scope :

Design and Development, Manufacture and Sales of Non-sterile, Non-active, Non-implantable Reusable Surgical Instruments and related Accessories

Technical Areas (Ref. IAF MD9):
 A.1.1

Re-certification : 16/09/2024
 Valid until : 17/09/2027
(Certificate validity is conditionalised by the fact that the previous certification cycle of IAF Accredited Certificate)

Authorised Signatory

The authenticity of this document may be confirmed by visiting certification@zenithquality.com. Lack of fulfillment of conditions as set out in the Certification Agreement may render this certificate invalid. Any alteration, forgery or falsification of the content or appearance of this document is prohibited and offenders may be prosecuted to the fullest extent of the law.

Accreditation Body: National Accreditation Board for Certification Bodies (NABCB)
 Quality Control of India (QCI), 2nd Floor, Institute of Engineers Building, Bharath Shree Park, Marol
 New Colony - 110 025, India
www.nabcb.co.in/qci/

Issuing Authority: Zenith Quality Assessors Pvt. Ltd.
 Office No. 12A, 4th Floor, Institute of Engineers Building, Bharath Shree Park, Marol
 New Colony - 110 025, India
 Email: info@zenithquality.com
www.zenithquality.com

CMMS/REG/2016-0000/10000 018 23

**Damaged Screw
 Extractor Set**

New!



Screw retriever



OR00329

Padmasree Dr.S.M.Hardikar
 Facilitates the removal of
 damaged screws. Suitable
 for 3.5 & 4.5 screws

**Locking Driver
 Left**



4.5 3.5 2.7

Carbide Bits



5 3.5 2.5 2

Hollow - mill



2.7 - OR00332
 3.5 - OR00334
 4.5 - OR00335
 6.5 - OR00336

**Dr. Mani
 JJ Bone Clamp**



OR00347 Forearm
 OR00349 Humerous
 OR01347 Paediatric

Designed by Prof. N.J. Mani
 A great time saver in Fore
 arm / Humerus reduction.
 Plate inserting facility without
 disturbing reduction

**Arjuna K Wire
 Cutter**



OR00377

K wire Cutter with tungsten
 Carbide Tip for K-wire
 up to 2.5 mm

**Sachin K Wire
 Cutter**



New!

OR00378

K wire Cutter with tungsten
 Carbide Tip for K-wire
 up to 1.5 mm

**Broken Interlocking
 Nail extractor**



OR00325

Unique thread design useful
 for Brocken nails. Designed
 to be used in a counter
 Clock wise direction

JJ Interlocking Rasp



OR00343
 Humerous

OR00342

Tiilbea & Femour

Dr.N.J.Karne,Pune
 The Cannulated Rasp to widen
 the interlocking entry point
 from 8mm to 13mm

**JJ Implant Extractor
 with Slap Hammer**



OR00338 Large
 OR00339 Small

The tungsten tipped extractor
 for square nails, interlocking
 nails,K-nails,K-wires,Rush nails
 knowled pins,AM pins, Screws etc.

**Wire Twister with
 Self Cutting**



OR00345

Wire Tensioner with automatic
 cutting. An innovating instrument
 for cerclage wiring.

**Universal Interlocking
 Nail Extractor**



9 Cannulated Adaptors

OR00330-03

Extractor for any interlocking
 nails of any make

Small bone Compressor/Distractors

New!



OR00412R-00
Radiolucent



With Fixed & Rotated Arm Alm- JJ



Fixed Arm- OR00412-21
Rotated Arm- OR00410-21



OR00411



OR00412

Dial Mechanism helps allow precise control of inserted wires in small bone surgery

Hintermann Distractor

Closed Arm
OR56438-1



Available in two sizes

To distract small joints in a linear direction in foot, hand and spine surgery



Outspreader Arm

OR56438-1OS



New!

To distract small joints in Episurf surgery

Calcaneal Spreader



RT418-01



Used for lateral column lengthening of the calcaneus

MIS Drill Guide Wire Passer



OR 8315-00

Very Useful for Wrist, Elbow, Foot-Ankle Reconstruction

Mippo Wire Passer



MP00001
For large bones



The semi sharp awl to make path around the Bone for passing mippo wire passer

Incision Aligner


New!



OR 1330
Aligning the instrument with Incision by Extending one hook outward. The sliding end locks in place when tensioned.

Bone Mill

Table mounted with crank handle



OR 6001-00
For grinding bone chips used for bone graft procedures

Nagpur Retractor



4 mm pin

OR4009

New!

Idea by Dr. Anil Golhar
 For verticle retraction conventional charny retractor is used over acetabulum & greater trochanter to spread in other direction this retractor can be used.

**Steinmann Pin Introducer
 Universal Chuck w/o Key
 Fully Stainless Steel**



OR 5507-14

14 cm/5 1/2" Length
 0.7 mm - 6 mm pins

**JJ Raindrop K-wire
 Holder**



OR00410

A Chuckless K-wire holder
 for 0.5 to 2.5mm K-wires

Gelpi-Blade Retractor



OR40080

Designed to help expose
 difficult to visualize areas at
 the end of incisions

Plate Holding Forceps



4.5
 OR6045-20

3.5
 OR6035-20

2.7
 OR6027-16

**Plate Holding Forceps with
 Screw Guide**



4.5
 OR60450

3.5
 OR60350

New!

Strayer Retractor



New!

9.25"

OR1869

Blade Length: 8.9cm

Blade Width: 1.5cm

To retract the soleus muscle for
 exposure of gastrocnemius fascia

Charney Retractor



OR4006-10



OR4007-75

OR4007-50



OR4006-75


OR4006-60

With two pairs of blade

Offset Caliper



New!



Offset Measurement



Femoral head measurement


OR 7334-30

Using a caliper to restore the centre
 of the femoral head during THR
 The device, which is adjustable for
 both leg length and offset, connects
 to this pin.

Double Barrell Drill sleeve for Distal Femur

Save Radiation Sleeve

New!



OR1110-01

Idea by
 Dr. Chandak, Nagpur

After fixing the first distal screw, You can fix the second screw without any "C" arm Exposure

PFN Entry Point Selector

New!



OR1050-01

Soft Tissue Bone Clamp

New!



10" - OR5047-25
 8" - OR5047-20

It is used for reduction of bone fracture with radiolucent soft tissue attachment.

Micro Fracture Drill Sleeve

New!



OR1106-100

OR1102-25

45° Angled

Double Barrell Femoral Guide Sleeve

New!



OR1110-02

Idea By : Dr. Adithya Bothra
 Femoral Guide pin insertion,
 pelvic and acetabulum fixation,
 reverse diamond fixation.

RadioLucent Sleeve



Custom design

Ulna Nerve Protecting Sleeve

New!



OR1108-25

Idea By : Dr. Shah Miraj
 This k-Wire sleeve to protect
 ulna nerve while doing k-wire
 fixation for supra condyle
 humerus in children

Pediatric Lateral Condyle Clamp

New!



OR5042-15

'S' Wire Bender



OR00375

Bending on 120°

K-Wire Bender/Cutter



OR01801

Can bend and cut
 K-Wires 1 to 1.6 mm



Bone Clamp with Drill Guide



HF502217

Chain Retractor

New!



Tendon Clamp with Slit

New!



- 1mm SHF 1.0
- 1.5mm SHF 1.5
- 2mm SHF 2.0
- 2.5mm SHF 2.5
- 3mm SHF 3.0
- 3.5mm SHF 3.5
- 4mm SHF 4.0
- 5mm SHF 5.0
- 6mm SHF 6.0
- 7mm SHF 7.0

For Cutting Tendon and Nerve Straightly

Scaphoid Lunate Clamps

New!



Reposition Forceps with K Wire Holder For Percutaneous Application

Reposition Forceps For Reduction of Hard to Reach Fractures

K Wire Inserter/Extractor

New!



For 1mm to 4.8mm Pins
 Adjust the Screw to Change the Size of K Wire

Guide Wire Holder



Quick holding and release of a Guidewire passing through the IM Canal. Also holding 1mm to 4mm K Wires

Asymmetric Webers Clamp

Straight, Curved, Reduction Forceps

New!



Make a 2.5mm hole in the shaft to accommodate the straight tip of the forceps curved tip will engage tip of Olecranon to allow placement of the plate while holding reduction. Also used medial malleolus fractures.

Femoral Tibial Spreader



Tibia Punch

New!



Designed to impact cancellous bone into the subchondral weight bearing region of the tibia

JJ Closed Reduction Forceps



OR00326

Idea-Dr. Harikumar, Kannur
 A radiolucent forceps for Closed reduction

Bone Awl 5"



Humrus Str. OR13614-02
 Humrus Cvd. OR13614-03
 Fore Arm Str. OR13614-00
 Fore Arm Cvd. OR13614-01
 Finger Str. OR13614-04
 Straight & Curved

Blade Retractors



15x5mm HF5215-05
 20x5mm HF5220-05 RT163
 15x8mm HF5215-08
 20x8mm HF5220-08
 Mini Langenbeck 5" Mini Sen muller 5"
 Available in any dimensions

Sharp Hook



OR00434

To Check Fracture and remove Ingrown tissue

Bone Rasp



OR6954

Bone Scoop



2mm OR723-02
 4mm OR723-04

Caroll Osteotome



2mm HF5291-02
 4mm HF5291-04
 6mm HF5291-06
 8mm HF5291-08
 9mm HF5291-09
 10mm HF5291-10
 12mm HF5291-12
 14mm HF5291-14
 16mm HF5291-16

Bone Graft Impactor



2mm OR723-02
 4mm OR723-04

Available in round and Rectangular shapes

Periosteal Elevator



3mmStr. HF7008-03
 5mmStr. HF7008-05
 3mmCvd. HF7018-03
 5mmCvd. HF7018-05

Straight & Curved
 3mm & 5mm

Tendon Aproximator



OR00408

Tendon Stripper 12" & 6"



HF5812-08
 HF5806-06

Available in two sizes

Tendon Seizing Forceps



HF5814-02

HF5810-15

Delicate

Tendon Interlacing Forceps



Available in
 4"
 5"
 6"

Straight & Curved

Tendon Passer



HF5805-25

Adson 1x1




GN2781

New!

Idea by : Dr. Sudheer Warriar

JJ Collinear Reduction Clamp

New! With New ratchet locking handle



OR 1-500
 A Minimally invasive reduction



Surgical experience by
 Dr. Sreehari (MIMS Kannur)

Double - Ended Hohman bone leavers



Cannulated Bone Awl



OR13112-00 OR13112

Suction tube and cautery holder



OR54330

JJ Open Reduction Forceps



OR00352

It allows accurate insertion of pins for internal fixation.

Distal Femur Distractor



OR4220-10

New!



OR4220-00

To distract distal femur from the proximal tibia.

Curette Osteotome



OR6908

Easy removal of osteophytes and cartilage without having to switch tools

Femoral Cement Osteotome



OR6907

To help trim away cement from around a femoral knee component

Charnley Ring Curette



T- handle
 OR 8430-37
 Double sided

37 cm/15" Length
 Serrated cutting edge

Intramedullary Irrigation Cannula

OR 1289-01



This cannula is very usefull to give a through wash in side the intramedullary canal

Osteotomy Impactor



20mm
 OR1850

New!

To impact hinge area without breaking

Nose Pliers



Available in Narrow nose, flat nose, long jaw needle & dolphin nose

Inter Changeable instruments



OR00317
 Osteotomes, gouges impactors, curettes, Rasps with quick coupling handle

Cannulated Wire passers



Wire passers for patella forearm humerus Femour and Tibia

Wire - Passers



Used for Forearm, Humerus Femour & Tibiya

Guide Wire Sterilization Container



OR00315
 Idea-Dr. Francis Roy
 Flexible, Autoclavable available in 1mtr. and 1.25 mtr.

Syndesmosis Clamp



Idea By Dr. Wolfson, USA

Condyle Clamp



Large OR721-14
 Small OR721-09
 Design by-Dr. Vijay Bose, Chennai

Condyle Clamp with Sleeve



OR722-00

JJ Bone Curette with Suction



OR0043-04

Pattela Holding Forceps



OR00363 OR00361 OR03610

Available with Ratchet lock Speed lock and Smart lock

Two pointed reduction Clamps



HF5042-23 HF5042-19 HF5042-14

Designed for Medical Malleolar / Bone Fragment

Bone Graft Holding Forceps



OR18502

Idea-Dr. Rajasekharan, Koval
 Bone graft can be grasped & placed

JJ Stem Extractor

New!



OR6909-25N

Additional Attachments
 Loop attachment
 Hook attachment

The instrument is used to extract stem during hip revision surgery.

Cement Grasping Forceps



OR6903-25

OR6901-35

Cement Removal Instruments



Overall Length 13"

Femoral Head Bone Plug Cutter




OR00409

Taking bone from femoral head to plug femoral canal before cementing

Flexible Osteotomes

7" OR3807-20 OR3807-15 OR3807-10 OR3807-05
 9" OR3809-20 OR3809-15 OR3809-10 OR3809-05



20mm

15mm

10mm

5mm


Sharp, Flexible blades are well suited for loosening implants from cement or bony in growth fixation

Ball Reamers



The extra long Ball Reamers may be used to ream cement layers as well as cement islands in the distal femoral cavity

Micro Forceps



6" Round body plain & Toothed

6" Flat body plain & Toothed

Castro viejo Needle Holder Scissors




STR MI610 MI520
 CVD MI611 MI521
 6" Straight & Curved

Osteophytes Osteotome



OR7619-18

Mippo Tunnelers



20cm

25cm

30cm

MP00020 MP00025 MP00030

Wire Holder



TW203



10mm jaw
 MI00419



20mm jaw
 MI00326



MI405

Hand Reamers
Cannulated Reamer



Rigid Reamers



Wire Passer/Retrieval Forceps



Capsule Scissors



**Anteversion
 Osteotome**



Pelvic Impactor




**Femoral Head
 Extractor**



**Sub Trochanteric
 Reduction Clamp**



Pelvic Reduction Forceps



Kolbel Self Retaining Glenoid Retractor with Center Blade



SM 1063-20

Center blade
 Depth 6.4cm/2.5"

SM 1063-21

Blades - pair
 36mm x 36 mm

SM 1063-36

Blades - pair
 36mm x 53 mm

SM 1063-53

Subacromion Spreader



SM 6552-19

SM 6551-19

Used for massive tears and repair of a supraspinatus tear

Femoral Distractors



OR4005-00

RadioGraphic Ruler



OR 5008-00

Multi purpose pliers
 With Strike Plate



Big OR00802
 Small OR00801 OR00802-00

Designed for bending and cutting of K-wires and plates. Available in two sizes

Plate Bending Plier



OR804

Distal Screw Aiming Device



OR04110

A Radiolucent device for distal locking aiming and drilling.

Plaster Cutting Guard

New!



Child Size
 OR1520-01

OR1520
 Designed to help protect the patient's skin when removing a plaster

JJ Ortho Suction Tube with Cleaning Brush



OR4230

Helps eliminate plugging due to bone, cement fragments blood clots etc...

Fracture Aligner



OR00327

Used to align the bone ends of a femur fracture.

New!



OR00328


For lesser Trochanter

Cement Removal



OR5245 OR5212 OR5211
 Double ended instrument designed for cement removal during arthroplasty surgery

Knee Self Retaining Retractors



Curved Hohmann, Narrow
OR 66500

Curved Hohmann, Wide
OR 66501

PCL Retractor
OR 3525

Patella Retractor
OR 66502

Collateral Ligament Retractor
Single Tip OR 66301

Collateral Ligament Retractor
Double Tip OR 66300

Z Knee Retractor



OR3111-11

Helps create better access to the articulating surfaces

Collateral Ligament Retractor



OR7661-10

Helps protect the lateral collateral ligament while exposing the proximal tibia

45° Knee Retractors



OR3108-01

OR3108-00

OR7742-20

Designed for general use around the knee

S" Total Knee Retractor



OR3111-10

Helps to protect the collateral ligaments and popliteal structures while providing excellent visualization within the knee joint

Honeycomb Drill Guide

OR221-00



New!

Available in Various Depths
100mm, 125mm & 145mm

Place the steinmann pin through the centre hole. If a slight adjustment is necessary, choose one of the alternate holes and slide the tissue protector off of the pin and back on to the new pin location with central hole aligned


Posterior Capsule and Sciatic Nerve Protection Retractors



Right : OR6115R
Left : OR6115L

Sits on the outer rim of the posterior/ inferior ischium to retract the posterior capsule for acetabular exposure and help to protect the sciatic nerve.

HIP Bone Levers



OR3108-45

OR3108-35

OR7780-20

OR7750-14

OR3108-24

OR3010-25

OR3112-40

OR3108-25

OR30285-34

OR30285-30

MIH Retractors



Overall Length 330mm
 Blade Width 25mm
 OR 2030



Overall Length 330mm
 Blade Width 32mm
 OR 2031



Overall Length 330mm
 Blade Width 40mm
 OR 2032



Deep Hohman
 Overall Length 330mm
 Blade Width 32mm
 OR 2033



Blunt Hohman
 Overall Length 330mm
 Blade Width 32mm
 OR 2034



Muller Hohman
 Overall Length 350mm
 Blade Width 25mm
 OR 2035

HIP Retractors



40° Angled
 Overall Length 280mm
 Blade Width 18mm
 OR 2451-28



Double Angled
 Overall Length 340mm
 Blade Width 41mm
 OR 2451-34



Blunt Retractor
 Overall Length 320mm
 Blade Width 24mm
 OR 2451-32



Blunt, Double Angled
 Overall Length 290mm
 Blade Width 24mm
 OR 2451-29



Bone lever After Ganz
 Overall Length 270mm
 Blade Width 24mm
 OR 2451-27



Pubis Retractor After Ganz
 Overall Length 260mm
 Blade Width 22mm
 OR 2451-26



Double Bend Hohmann
 Overall Length 280mm
 Blade Width 15mm
 OR 2451-28



Double Bend Hohmann
 Overall Length 310mm
 Blade Width 22mm
 OR 2451-31



Double Bend, Blunt Hohmann
 Overall Length 300mm
 Blade Width 20mm
 OR 2451-30



Bone Lever Amin Wide
 Overall Length 320mm
 Blade Width 32mm
 OR 2451-32



Bone Lever 65° V-Shaped Spike
 Overall Length 320mm
 Blade Width 30mm
 OR 2451-70



Trochant Bone Lever V-Shaped Long Spike
 Overall Length 330mm
 Blade Width 30mm
 OR 2451-33

HIP Retractors



Double Bend 2 Prongs Hohmann
 Overall Length 300mm
 Blade Width 17mm
 OR 2451-76



Double Bend, Wide, Blunt Hohmann
 Overall Length 310mm
 Blade Width 32mm, Tip Width 15mm
 OR 2451-78



Blunt Hohmann
 Overall Length 280mm
 Blade Width 32mm, Tip Width 15mm
 OR 2450-50



Neck Elevator, Muller
 Overall Length 380mm
 Blade Width 26mm
 OR 2451-90



90° Angled Deep Retractor
 Overall Length 210mm
 Blade Width 12mm
 OR 2452-21



Double Bend Hohmann Retractor
 Overall Length 303mm
 Blade Width 15mm
 OR 6210-02



Pelvis Retractor
 Overall Length 310mm
 Blade Width 22mm
 OR 1310-70



Cobra Retractor
 Overall Length 340mm
 Blade Width 36mm
 OR 1310-60



Cobra Retractor
 Overall Length 340mm
 Blade Width 30mm
 OR 1310-00



Bent Hohmann Retractor
 Overall Length 310mm
 Blade Width 22mm, Sharp Tip
 OR 1310-40



Bent Hohmann Retractor
 Overall Length 310mm
 Blade Width 22mm, Blunt Tip
 OR 1310-50



Bent Hohmann Retractor
 Overall Length 310mm
 Blade Width 22mm, Stub Tip
 OR 1310-10

Bent Hohmann Retractor 90°
 Overall Length 250mm
 Blade Width 19mm
 OR 3195-19



Bent Hohmann Retractor 90°
 Overall Length 185mm
 Blade Width 19mm
 OR 3108-19



AVA Retractor
 Overall Length 285mm
 Blade Width 22mm, Blunt Tip
 OR 4540-01



AVA Retractor
 Overall Length 285mm
 Blade Width 22mm
 OR 4540-02



Femoral Neck Elevator
 Overall Length 388mm
 Blade Width 32mm, Cup Model
 OR 1309-40



Femoral Neck Elevator
 Overall Length 325mm
 Blade Width 25mm
 OR 1309-50



Femoral Neck Elevator
 Overall Length 332mm
 Blade Width 25mm
 OR 1309-70

HIP Retractors with Grip Handle



HIP Skid
 Overall Length 300mm
 Blade Width 29mm
 OR 1309-20



Cobra Retractor
 Overall Length 340mm
 Blade Width 30mm
 OR 1310-10

New!



Pelvis Retractor
 Overall Length 310mm
 Blade Width 22mm
 OR 1310-71



Double Bend Hohmann Retractor
 Overall Length 390mm
 Blade Width 15mm
 OR 1308-51



Bent Hohmann Retractor 90°
 Overall Length 340mm
 Blade Width 19mm
 OR 1308-71



AVA Retractor
 Overall Length 370mm
 Blade Width 22mm, Blunt Tip
 OR 1309-11

Radiolucent Bone Levers

New!




Cobra Retractor
 OR 1310-30R



Cobra Retractor
 OR 1310-35R



Shallow Bone Lever
 OR 7755-25R



Succesfull Feedback Reference

- Dr. Babu Joseph**
(Renai Medicity, Kochi, 9847165295)
- Dr. Karthik Shukla**
(Ahemadabad, 9825296189)
- Dr. Dr. Rajeev Raj Choudhari**
(Surat, 9824057579)
- Dr. Dr. Swarnendu Samanta**
(Kolkatha, 9830083031)

Radiolucent Bone Levers

New!



Modified Hohmann
 OR 4540-20R



22mm Standard Hohmann
 OR 3018-22R



30mm Standard Hohmann
 OR 3018-30R



Humeral Retractor
 SM 0592-21R



Browne Deltoid Retractor
 SM 0500-21R

Radiolucent Retractors

New!

30 x 10



OR2530-10R

30 x 30




OR2530-30R

50 x 30



OR2550-30R


CASPAR MICRO LUMBAR DISCECTOMY RETRACTOR




Micro Discectomy Retractor PS5200-00



Caspar Vertebral Body Distraction Set



Super -Slide Laminectomy Retractor




Laminectomy Retractor



Finochitto retractor with hinged arms and detachable blades.

Cervical Retractor set



Cobbs Periosteal Elevator



Bone Marrow Aspiration Needles



Bone Biopsy Needles with Trocar and Ejector



Bone Graft Funnel



Kerrison Laminectomy Rongeurs



Bayonet Cervical Curette



Foramenotomy Punch



Ferris-smith Kerrison Punch



Black Plasma Coated, thin Foot Plate, detachable with ejector

Long Anterior Kerrison Punch



Intervertebral Disc Rongeur



Kerrison Punch Down



Intervertebral Disc Rongeur Black



JJ Bone Clamps with Drill Guide



- OR1890-02 9.25"
Large- 10mm Drill Guide
- OR1890-01 6"
Small- 8mm Drill Guide
- OR1890-00 4"
Extra Small- 6mm Drill Guide

Designed for the placement of bone clamp where it can best stabilize bone fragments. Drill guide allows screw placement through top of clamp.

JJ Modified Meniscal Clamp




- OR1883
Length: 7" (17.8 cm)
Jaw length: 1.5" (3.8 cm)

Designed to help facilitate grasping the posterior portion of meniscus. Improved bite when tension is placed on the meniscus and also used to help remove fat pad and suprapatellar bursa.




JJ Femoral Neck Cutting Guide



- OR4590
Overall Length: 8.25" (21 cm)
Shaft angle at end: 30°

Designed to help determine the femoral neck osteotomy location during the anterior approach for THA.


JJ Fracture Distractor



- OR1882
Overall length: 17.8 cm

Designed for providing opposing leverage to help bring the fibula (or any bone) back out to its proper length after it has been shortened by a fracture.

JJ Hip Length Gauge



- OR1327-00

Helps to determine the leg length and hip offset in total hip arthroplasty, including minimally invasive techniques.

JJ Malleable Bone Tamps



- OR5296
Overall length: 37 cm
Ø 13 mm
- OR5296-01
Overall length: 25 cm
Ø 10 mm

The large tamp helps to elevate a depressed tibial plateau fracture. The small tamp helps to elevate a depressed tibial plateau and smaller tibial plateau fractures.

JJ Modified Pin Inserter



- OR3016
K-wires upto: 1.4 mm
Tube Length: 47 mm
Overall Length: 11.4 cm
Handle Length: 9.8 cm

Designed to aim and control the placement of flexible k-wires when they contact hard cortical bone, while helping to protect neurovascular structures from the spinning wire.

JJ Tibia Punch



- OR5140
Overall Length : 17.8 cm

Designed to impact cancellous bone in the subchondral weight bearing region of tibia. This helps to improve mechanical interlock in the cancellous bone/cement interface.

JJ Tibial Plateau Elevators



OR5298 OR5297
 Finish Starter
 Ø 10.4 mm Ø 4.7 mm
 Overall Length
 27.9 cm
 Idea By : Dr. Chandak,
 Nagpur

Designed to help with indirect reduction of a depressed tibial plateau fracture, and can be used with arthroscopic visualization and percutaneous fixation.

JJ Osteophyte Grasper



OR1776
 Overall Length: 11.20" (28.5 cm)
 Shaft Length: 7.28" (18.5 cm)
 Shaft Width : 7 mm
 Jaw Width at End: 4 mm
 Toothed Jaw Length: 14 mm

Useful in helping to remove posterior osteophytes in knee surgery, and helps to remove the labrum and soft tissues in anterior total hip surgery.

JJ Patella Double Scalpel Handle



GN138-06
 Length: 5.75" (14.6 cm)
 #3 handle size

The offset of blade is 10mm. Harvests tendon graft from the patella and tibial tubercle including the patellar tendon.

JJ Pin Clamp



OR1760-01
 Overall Length: 6"
 1.8 mm cannula

Designed to allow accurate insertion of pins for internal fixation.

JJ Ruler With Right Angle Handle



OR1450
 Handle Length: 4.25" (10.8 cm)
 Ruler Dimensions: 2.5 x .5"

90° angle handle. Designed to measure the femoral head/neck length.

JJ Ruler with 45° Angle Handle



OR1430
 Handle Length: 5" (12.7 cm)
 Ruler Dimensions: 2.5" x .5"

45° angle handle. Designed to measure distances in small deep incisions.

JJ Vernier Caliper



OR7338-00
 Length 22.5 cm
 Measuring Range 15 cm

This instrument is used for measuring the disc space and diameter of the femoral head.

JJ Small Bone Tamp with Oblique K-wire Hole



OR5294
 Wire hole 1.1 mm k-wire
 Overall length: 19.1 cm
 Shaft diameter: 6.3 mm
 End diameter: 2.5 mm

OR5294-01
 Wire hole 1.6 mm k-wire
 Overall length : 19.1 cm
 Shaft diameter : 6.3 mm
 End diameter : 2.5 mm

This design helps in the concurrent reduction of a fracture and placement of a wire into the fracture site. It is particularly helpful when the surgical exposure is small and tight, the fracture fragments are small, and the reduction is demanding.

JJ Radiolucent Retractors



OR22550-R2K
 length 21.5cm
 Ø3.2 mm K-wires

OR25230-R
 length 24 cm

OR25210-R1K
 length 20cm
 Ø 3.2 mm K wire

OR22550-R2K: Designed to retract the iliac vessels.

OR25230-R: Acetabular Rim Retractor. To be placed near ilio-pubic eminence.

OR25210-R1K: To be placed near pubic tubercle.

JJ 90° Bone Hook



OR5940-S Sharp Tip
 OR5940-B Blunt Tip

Designed to ergonomically help the surgical assistant elevate the proximal femur during TKA.

JJ Femur/ Tibia Fracture Distractor



OR1809
 Overall Length : 10.5" (26.7cm)
 Overall Width : 7.25" (18.4cm)
 For Pins upto : 0.25" (6.4mm)

JJ Femur/ Tibia Fracture Distractor is used with self centering bone clamps for diaphyseal fractures and 6mm Schanz pins to distract intra-articular fractures for reduction and fixation.

JJ Anterior Hip Referencing Rod Assembly



OR2674-00 - Full assembly
 Overall Length: 27.75"
 Rod Diameter: 0.25"

OR2674-A - Top assembly
 Overall Length: 16.75"
 Rod Diameter: 0.25"

OR2674-B - Bottom assembly
 Overall Length: 10.5"
 Rod Diameter: 0.25"

It helps to determine implant fit, position, alignment and recreation of leg length and offset using the contralateral hip for reference.

JJ Goniometer with Ruler



OR2029

JJ Goniometer is used to measure an angle or to permit rotation of a joint to a definite position. Ruler attachment helps in precise measurement along the area.

JJ Bayonet Mani Clamp



Radius
 Right : ORB00347R
 Left : ORB00347L
 Humerus
 Right : ORB00349R
 Left : ORB00349L

An updated bayonet design of JJ Mani clamp.

JJ Radiolucent Awl



OR8030-01

JJ Radiolucent Owl is used to locate the holes in interlocking nails.

JJ Amputation Shield



OR8898

This instrument is designed to use as a shield for protection during amputation procedures.

JJ Adjustable Knee Positioner

OR2770-00



Successfull Feedback
 By
 Dr. Hemant Maurya
 Shikha Hospital, Allahabad



Height:16" Height:14" Height:11"

Available in Fixed Model

JJ Adjustable Knee Positioner with Leg Support

OR2770-LS



New!

JJ Adjustable Knee Positioner is used for proccures around the knee. The adjustable design helps in positioning at different heights and angles.

JJ Adjustable Knee Positioner with leg support for additional support for the foot.

JJ Acetabular Knife



OR11122
 Length : 10"

JJ Universal Tibia/ patella clamp

OR3685
 Length : 10"



The design helps in removal of tibia wedge and also helps in everting patella.

Malanski Retractor



OR3110-30
 Overall Length: 9"

Whelan Retractor



OR7116
 Blade Depth: 6"
 Overall Length: 14.5"

JJ Canal Opener



OR3006-01
 Overall length : 60cm

JJ PFN Aimer

OR5008



Ø3 mm Holes

Designed by
 Dr. Harjot Singh Gurudatta

Contributions by
 Dr. RM Chandak &
 Dr. Abhishek Tiwari

Ø3 mm Holes

This device is a PFN Aimer which helps in determining the cranial and caudal placement of two screws in two screw PFN system. This is a radiolucent device and will avoid the multiple k wire insertion and multiple adjustments while insertion of a PFN.

JJ Ratchet Bone Clamp

OR6180-23R

Overall length : 23.5cm
 Jaw Opens to : 8 cm



JJ Ratchet Bone Clamp is designed for quick tightening and release around the bone.

T Handle Ruler



OR443-040
 20mm

New!

To make measurements inside the surgical site to assist in implanting a prosthetic femoral hip into the intramedullary canal of the natural femur.

JJ Chain Retractor



RT6828-00 Complete Set

Four Retracting Blades
 Blade Holder
 Mounting Frame

JJ Chain Retractors is a self retaining table mounted retractor used for General and Gynaecology procedures. The set includes four retracting blades and a blade holder along with chains which helps in effective retracting at different heights.

JJ Thomson Retractor

TRT6969



JJ Thomson Retractor is designed to provide access and exposure for various procedures such as all Spine Surgeries, Transplant Surgeries, Abdominal General Surgery and Abdominal Vascular Surgery.

Last Screw Retractor


New!



4" depth variant with 30° for pelvis



OR1764-30



OR1766
 Length : 8.5"

It is used to widen the incision to insert the last screw.

Cobbs Elevator with Suction



OR7671-17
 Length : 13"

New!



It is used to provide better exposure during spine, hip and trauma cases with a suction provision.

Dissectors with Suction

Penfield Dissector With Suction
 Length : 7"
 PS9406



New!

Dura Retractor with Suction Frazier
 PS5013-19
 Length : 8"



Dissectors used for elevating with a suction provision.

JJ Weitlaner Retractors with Radiolucent Tip

New!



RT4417-33R
 Length : 7"

Weitlaner Retractor with Radiolucent tips used for retraction of tissues and muscles.

JJ Gelpi Retractors with Radiolucent Tip

Gelpi Retractor

New!



OR4030-16R
 Length : 7"

Gelpi Retractor with Radiolucent tips used for retraction of tissues and muscles.

Gelpi Retractor with Flat Tip

New!



OR4400-24R
 Length : 10"

Gelpi with Radiolucent Flat tips used to provide better holding grip during retraction.

Titanium Needle Holder

New!



Length : 8"
 NH153T

Needle Holder made of Grade 5 Light weighted Titanium.

Titanium Weitlaner Retractor

New!



Length : 5"
 RT4335-23T

Weitlaner Retractor made of Grade 5 Light weighted Titanium.

Titanium Vascular Clamps

New!



Vascular Clamps made of Grade 5 Light weighted Titanium.

JJ Interchangeable Impactor Set



- Round Head 09mm
- Round Head 09mm
- Square Head 10 x 10mm
- Rectangle Head 12x 08mm
- Oval Head 13 x 08mm

Length : 9"

OR11122

New!

The Special JJ Interchangeable impactors Set with five different impactor heads.

JJ Humeral Head Holding Forceps Wire Guide



New!
Length : 9"

JJ Humeral head holding forceps with three 3.2mm wire hole on each side.

JJ Glen Graft

OR401-531



New!

In a Glenoid shoulder recon, when an offset or angled glenoid is discovered, the surgeon may need to place a graft to "level" the surface for the implant. The instrument we are creating is a graft holding jig. The angles of K wire placement are key as well as the centering jig for the implant tap piece.

JJ Radiolucent Finger Positioner



OR1134
Overall Length : 7"
Handle Width : 4"
Blade Width : 40mm



New!

The fully radiolucent JJ Finger Positioner is used to provide the required positioning of finger during hand surgical procedures.

JJ Radiolucent Murphy Skid



OR0589-23R
Overall Length 11"

New!

The fully radiolucent JJ Murphy Skid used to retract muscle and tissue.

Flexible Reamer



OR69020

New!

The flexible reamer system is intended to clear the medullary canal of bone marrow and debris and to effectively enlarge the medullary canal of long bones for the acceptance of an intramedullary nail.

MCL Retractor



Length : 8.625"
OR7109

New!

It is used for retracting the femur and protect the MCL.

JJ Carpel Tunnelor Guide

HF5888-03

New!



JJ Carpel Tunnelor with tissue protector and knife.

Versatile Small Bone Spreader



OR1829

Overall Length : 5"

Blade Depth : 28mm

Blade Tapers from 10mm to 6mm

New!

The small tapered blade helps in retracting small tight bone areas.

JJ Shoulder Retractor Frame and Blade Set

OR4696-00



New!

Blade Size 41mm x 55mm



Blade Size 43mm x 45mm



Blade Size 32mm x 35mm



228mm x 203mm

JJ Distal Radius Mini Frame and Blade Set

OR1578-00



New!

Blade Size 15mm x 15mm




Blade Size 15mm x 25mm



Blade Size 15mm x 35mm



Frame 70mm x 60mm



JJ Radiolucent Distal Radius Mini Frame and Blade Set

OR1578-00R



New!

Blade Size 15mm x 15mm



Blade Size 15mm x 25mm



Blade Size 15mm x 35mm



Frame 70mm x 60mm



JJ Needle Nosed Vise Grip Parallel Plier with Soft Lock



OR2008-25
 Length :9"

New!

JJ special needle nosed vise grip plier with parallel opening and soft lock.

Laterjet Instruments

Coracoid Drill Guide



New!

OR7000-17

The Coracoid Drill Guide is used to shape and size the coracoid and make two drill holes in the coracoid.

Laterjet Screw Insertion Guide



OR107-12
 6mm Offset

New!

It is used for laterjet screw insertion during glenoid surgery.

Shoulder Retractor Set *New!*

Dual Point Glenoid Retractor

SM861-104



10"
 20mm

Two Prong Forked Retractor

OR5842



10.5"
 20mm

Hohmann Retractor

OR3019-20



10.25"
 20mm

Wolfe Retractor

OR861-108



11"



6 Teeth

Bhattman Retractor

OR861-004



9"

Deltoid Retractor

SM1001-50



Length 10.5"
 Depth 2"

JJ Nathanson Retractor

RT7788



New!

Small

Medium

Large

Extra Large

It is used for retracting the liver during laparoscopic surgery.

JJ Broken Screw Removal Instruments



New!

OR1329-00

Extended broken screw removal set

Toronto Mallet



New!

OR3900-00A

Toronto Mallet made of delrin for soft impacting.

Bone Screw Ruler Gauge



New!

OR443-752

It is used to measure screw length at table top and also to measure k wire sizes.

K Wire Sizer



New!

KSW001

It is used to measure various k wire sizes from 0.7mm to 3.0mm.

Tibial Plateau Forceps




New!

OR5299

Tibial Plateau forceps is used for grasping tibial bone fragment.

JJ White Line Marker

Idea By
 Dr. Ayman Ahmed Hilaly
 Khoula Hospital, Oman



OR2029-00

New!

The modified design of goniometer which helps to intraoperatively measure the angle between white side line and posterior condylar line.

Femur Tibia Spreader with Soft Pad



New!

OR2020-00R

Femur tibia spreader with soft pad attachment is used to separate femur and tibia.

Scarf Bone Holding Forceps



New!

OR6150-17
 Length : 7"

It is used for holding small bone during surgeries.

Kolbel Retractor with Radiolucent Blades



New!

SM1062-23R
 Blade depth
 53mm

Kolbel Retractors with radiolucent blades are used to retract tissue during shoulder surgery.

Revision Component Remover



New!

Standard
 OR5120-01

Offset
 OR5120-02

The instrument is specially designed for the extraction of components during revision surgery.

Carpel Tunnelor

New!



Bayonet Shaped

HF5889-00

It is used as tunnelor for metacarpel.

Malleable Shoulder Retractor

New!



Length : 15.5"

Deep
22mm Width
OR4537-01

Shallow
22mm Width
OR4537-02

Wide
30mm Width
OR4537-03

Malleable Shoulder retractors are used to provide better exposure during surgery.

JJ Muller Compression Device

New!



OR03490

It achieves interfragmentary compression by tightening a tensioner that is temporarily anchored to the bone and the plate.

V Retractor

OR4851

New!



Overall Length : 9"

Depth 7cm

It is useful in providing deep access to the glenoid rim during shoulder surgery.

Elbow Retractor

New!

SM4697

Length: 7"



Designed to retract the anterior capsule in the elbow for better exposure.

Deltoid Elevator

OR5081

New!



It is used for lifting the deltoid for better exposure of the Subacromial Space.

Radiolucent Shoulder Slide

OR6879R

Length: 11"



New!

Radiolucent Shoulder Slide is used to protect the prosthetic bearing surfaces during shoulder arthroplasty.

New!

JJ Tension Retractor



OR1586

Width x Depth

30mm x 75mm

30mm x 50mm

30mm x 25mm



Length : 9.5"

JJ Tension retractor with different blades provides enhanced exposure during surgical procedures.

JJ Lonstar Retractor

New!



LS302-113

LS302-220

Hook attachments

Rake Hooks

Sharp Hooks

Blunt Hooks




Biocompatible, Reusable, Radiolucent & Light weight.

Instruments with Spring and Ring Handle

Forceps
SRF-13A



Scissors
SRS-13A



Idea By
Dr. Nabarun Saha
Berhampore

JJ Interchangeable Screwdriver Set

New!



ORD7017-10

- Quick Coupling Handle**
 - Straight Model
 - T handle Cannu. Model

Screw Driver Heads

- Hex 1.5
- Hex 2.0
- Hex 2.0 Cannu.
- Hex 2.5
- Hex 2.5 Cannu.
- Hex 3.5
- Hex 3.5 Cannu.
- Star T6
- Star T8
- Star T9
- Star T10
- Star T25
- Flat
- Phillips

Set with 2 Quick Coupling Handles
and 14 Screwdriver Heads

JJ Radiolucent Light Weight Cervical Retractor

New!
RT1579-01R



- Blades Depth**
 40mm
 50mm
 60mm
 70mm

JJ Parallel Drill Sleeve

Adjustable




Available in
65mm OR1021-65
75mm OR1021-75

New!

Parallel Drill Sleeve with
adjustable sleeve position.

JJ Double Ended Radiolucent Awl

With 4 Interchangeable Pins
OR8030-02



New!

- 2mm
- 3mm
- 4.5mm
- 5.5mm

JJ Radiolucent double ended awl with
four interchangeable pins is used to
locate the holes in interlocking nails.

JJ Fossa Finder with Handle

OR1000



New!

Used to locate the piriformis
fossa or greater trochanter for
guide wire insertion.

JJ Oval Retractor

Radiolucent



New!

HF1580-00R
Frame & 5 Pairs of Retractor Hooks
The light weight retractor set is
used for hand surgeries

**JJ Strong Extractor
Forceps with T Handle**

3.5
ORD00329-35T

4.5
ORD00329-45T



New!

Strong extractor forceps for
removal of screws.

JJ Chuck with Delrin T Handle

OR5505-014

With Key

New!



Chuck with key for holding k wires.

Bone Holding Forceps Ikkuta



New!

HF5055-20
Long Ratchet

It is used in Shoulder Surgery

JJ Femoral Component Trunion Clamp

OR6911-23




Bayonet
9"

New!

JJ Acetabular Screw Holding Forceps



New!



Angled Up

OR1490-25

JJ Reduction Clamp with Drill Sleeve



New!

8"
OR5042R-21S
9"
OR5042R-23S

Long Ratchet

JJ Light Weight Bone Levers with Handle

Superior Femoral Neck Elevator

OR1440-07



New!

Various bone levers available with this design.

JJ Bone Levers with Silicon Moulded Handle

Double bent two prong retractor with serrated tip



New!

Ergonomic Silicon Handle

Various bone levers available with this design.

Glenoid Exposure Spreader

Left- SM6551-20L

Right - SM6551-20R



New!

Used in shoulder surgery for exposure of glenoid.

JJ Bone & Plate Holding clamp with Swivel Head



New!

Swivel Head
3.5
OR6035-10
4.5
OR6045-10

Bone and Plate Holding forceps with tiltable Swivel head for better holding.

JJ Parallel Tibial Spreader with Soft Lock



Parallel Spread upto 25mm

New!

OR2025-14S

Used for Parallel Spreading of Femur Tibia

Femur Tibial Spreader with Silicon Pad



New!

OR2024-17SP

Femur tibia spreader with silicon pad attachment is used to separate femur and tibia.

Acetabular Lining Impactor



New!

PS6918-21

JJ Broach Extractor

OR6909-23

New!



With Slaphammer and Rod

It is used to extract the hip broach when the broach post is broken.

JJ Cup Liner (Shell) Extractor



New!

OR6909-46

It is used for extraction of cup liner or shell.

JJ Die-Punch Fracture Elevator

OR617-00

New!



Elevator for die-punch fracture.

JJ Wrist Drill Guide

HF1333

New!



Wrist drill guide for TFCC Ulnar Tunnel Repair.

Double Curved Periosteal Elevator

New!

OR7020D-04S



To harvest bone graft from bones

JJ Bone Graft Harvester

OR03336-08 8mm
 OR03336-10 10mm
 OR03336-12 12mm

New!



JJ Fracture Compression clamp

New!



On Bone



With Plate

8" OR5046-20

9" OR5046-23

Used for compression of bones with and without plates.

JJ Malleable Curette

Idea By Dr. Vikas Gupta

New!



Available in 10.5" of various cup sizes.


Instruments with Spring and Ring Handle

Forceps

SRF-13A



Idea By Dr. Nabarun Saha Berhampore



New!

Scissors

SRS-13A



New!

JJ THR Square Head Pin Retractor



OR1116-35

New!

It is used in THR surgery for retraction.

Carpel Tunnel Retractors

RT1730-15

RT4336-56CT

New!



New!

They are used in carpel tunnel surgery.

JJ 'One Stop Store' for all Surgical Instruments

Orthopaedic Surgery

International Instruments
 Designers & Manufacturers of Surgical Instruments

Orthopaedic Surgery

USFDA, CE & ISO REGD.

Kinfra Park P.O., Koratty, Thrissur District, Kerala, India 680309 | +91 480 2735959, 2735858, +91 94970 77555
 sales@myjonline.com | myjonline@yahoo.com
 www.myjonline.com
 INDIA+USA+OMAN+UAE

Spine Surgery

International Instruments
 Designers & Manufacturers of Surgical Instruments

Spine Surgery

USFDA, CE & ISO REGD.

Kinfra Park P.O., Koratty, Thrissur District, Kerala, India 680309 | +91 480 2735959, 2735858, +91 94970 77555
 sales@myjonline.com | myjonline@yahoo.com
 www.myjonline.com
 INDIA+USA+OMAN+UAE

Cardiovascular Surgery

International Instruments
 Designers & Manufacturers of Surgical Instruments

Cardiovascular Surgery

USFDA, CE & ISO REGD.

Kinfra Park P.O., Koratty, Thrissur District, Kerala, India 680309 | +91 480 2735959, 2735858, +91 94970 77555
 sales@myjonline.com | myjonline@yahoo.com
 www.myjonline.com
 INDIA+USA+OMAN+UAE

General Surgery

International Instruments
 Designers & Manufacturers of Surgical Instruments

General Instruments

Orthopaedic | Spine | Hand & Micro
 Plastic | Maxillofacial | Laparoscopy
 General | Gynec & Obstetrics | Uro
 E.N.T | Cardiovascular | Thoracic | Neuro

USFDA, CE & ISO REGD.

Kinfra Park P.O., Koratty, Thrissur District, Kerala, India 680309 | +91 480 2735959, 2735858, +91 94970 77555
 sales@myjonline.com | myjonline@yahoo.com
 www.myjonline.com
 INDIA+USA+OMAN+UAE

JJ 'One Stop Store' for all Surgical Instruments

Plastic Surgery

INTERNATIONAL INSTRUMENTS
 Designers & Manufacturers of Surgical Instruments

Plastic Surgery

JJ
 always innovative.

USFDA, CE & ISO REGD.

Neuro Surgery

JJ International Instruments
 Designers & Manufacturers of Surgical Instruments

Neuro Surgery

USFDA, CE & ISO REGD.

Kinfra Park P.O., Koratty, Thrissur District, Kerala, India 680309 | +91 480 2735959, 2735858, +91 94970 77555
 sales@myjonline.com | myjonline@yahoo.com
 www.myjonline.com
 INDIA • USA • OMAN • UAE

Laparoscopy Surgery

JJ International Instruments
 Designers & Manufacturers of Surgical Instruments

Laparoscopy Instruments

Orthopaedic | Spine | Hand & Micro
 Plastic | Maxillofacial | Laparoscopy
 General | Gynec & Obstetrics | Uro
 E.N.T | Cardiovascular | Thoracic | Neuro

USFDA, CE & ISO REGD.

Kinfra Park P.O., Koratty, Thrissur District, Kerala, India 680309 | +91 480 2735959, 2735858, +91 94970 77555
 sales@myjonline.com | myjonline@yahoo.com
 www.myjonline.com
 INDIA • USA • OMAN • UAE

Gynec Surgery

JJ International Instruments
 Designers & Manufacturers of Surgical Instruments

Gynec Surgery Instruments

Orthopaedic | Spine | Hand & Micro
 Plastic | Maxillofacial | Laparoscopy
 General | Gynec & Obstetrics | Uro
 E.N.T | Cardiovascular | Thoracic | Neuro

USFDA, CE & ISO REGD.

Kinfra Park P.O., Koratty, Thrissur District, Kerala, India 680309 | +91 480 2735959, 2735858, +91 94970 77555
 sales@myjonline.com | myjonline@yahoo.com
 www.myjonline.com
 INDIA • USA • OMAN • UAE



International Instruments

Designers & Manufacturers of Surgical Instruments

US FDA, CE & ISO REGD.

ABOUT US

JJ International Instruments - A leader in designing and manufacturing of high quality stainless steel surgical instruments in India. With a vast experience of 46 years in the industry, JJ provides unmatched service and unbeaten quality to the worldwide market. Corporate office and manufacturing unit with state of art facilities is located in south Indian state, Kerala - The renowned 'God's Own Country'.

CONTINUES SUCCESS...!

JJ is now worldwide since exhibited in AAOS 2013 with tremendous response from both domestic and international surgeons. JJ now participates several tradeshows in USA and more than 30 tradeshows all over India yearly.

" International Quality @ Indian Price "

CE Reg No. 2P191227.IIWT46
ISO Reg No. MD-QMS/91/R/N/2643
US FDA Reg No. 3010077069
CDSCO Reg No:
Intern-Thris-KL/M/MD/010840

INDIA • USA • OMAN • UAE

JJ design & manufacture complete Surgical Instruments for - General Orthopaedic, Spine, Hand & Micro Surgery, Gynec & Obstetrics, Uro, E.N.T, Neuro, Thoracic & Cardiovascular, Plastic, General surgery.

Features and Benefits

- Genuine 'JJ' Brand
- Custom made instruments
- Made out of German Stainless Steel- SS 410 & SS 420
- Follows International Passivation standard - ASTM 967
- International hardness of product - 45 HRC & 50 HRC
- Universal Specifications & Quality
- 5 years replacement warranty against manufacturing defect
- Absolute Competitive price
- JJ welcomes innovative ideas and designs
- Gift packs on attractive terms

JJ's repairing & refurbishing program:

- Repairing & refurbishing of surgical instruments of any brand
- Reduce purchase budget drastically
- Annual Maintenance Contract (AMC)
- Round the clock service center
- Replace damaged instrument with new until rectification

JJ's Mobile Trade Expo

- Introducing mobile sales and service all over India
- Specially designed instruments on display
- Easy to select, exchange and experience
- Get technical advice from the experts
- Contact us for mobile trade expo schedules

www.myjjonline.com

sales@myjjonline.com, myjjonline@yahoo.com

